


Product Name: Recombinant Human NOVA1 Protein
Catalog #: PEH8281



Summary

Name	NOVA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P51513
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTHTGVPIDL DPPDSRKRPL EAPPEAGSTK RTNTGEDGQY FLKVLIPSYA AGSIIGKGGQ TIVQLQKETG ATIKLSKSKD FYPGTTERVC LIQGTVEALN AVHGFIAEKI REMPQNVAKT EPVSILQPQT TVNPDRIKQT LPSSPTTTKS SPSDPMTTSR ANQVKIIVPN STAGLIIGKG GATVKAVMEQ SGAWVQLSQK PDGINLQERV VTVSGEPEQN RKAVELIIQK IQEDPQSGSC LNISYANVTG PVANSNPTGS PYANTAEVLP

image



Background

Alternative Names	RNA-binding protein Nova-1 (Neuro-oncological ventral antigen 1) (Onconeural ventral antigen 1) (Paraneoplastic Ri antigen) (Ventral neuron-specific protein 1)
--------------------------	---

Product Name: Recombinant Human NOVA1 Protein
Catalog #: PEH8281



Note

For Research Use Only , Not for Diagnostic Use