


**Product Name: Recombinant Human NOVA1 Protein**  
**Catalog #: PEH8282**



## Summary

<b>Name</b>	NOVA1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P51513
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	29 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RKAVELIIQK IQEDPQSGSC LNISYANVTG PVANSNPTGS PYANTA EVLP TAAAAAGLLG HANLAGVAAF PAVLSGFTGN DLVAITSALN TLASYGYNLN TLGLGLSQAA ATGALAAAAA SANPAAAAAN LLATYASEAS ASGSTAGGTA GTFALGSLAA ATAATNGYFG AASPLAASAI LGTEKSTDGS KDVVEIAVPE NLVGAILGKG GKTLVEYQEL TGARIQISKK GEFVPGTRNR KVTITGTPAA TQAAQYLITQ RITYEQGVRA

## image



## Background

<b>Alternative Names</b>	RNA-binding protein Nova-1 (Neuro-oncological ventral antigen 1) (Onconeural ventral antigen 1) (Paraneoplastic Ri antigen) (Ventral neuron-specific protein 1)
--------------------------	---

**Product Name: Recombinant Human NOVA1 Protein**  
**Catalog #: PEH8282**

---



**Note**

For Research Use Only , Not for Diagnostic Use