


Product Name: Recombinant Human SH3BP2 Protein
Catalog #: PEH8289



Summary

Name	3BP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78314
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PMKAIGAQNL LTMPGGVAKA GYLHKKGGTQ LQLLKWPLRF VIIHKRCVYY FKSSTSASPQ GAFSLSGYNR VMRAAEETTS NNVFPFKIIH ISKKHRTWFF SASSEERKS WMALLRREIG HFHEKKDLPL DTSDSSSDTD SFYGAVERPV DISLSPYPTD NEDYEHDDDED DSYLEPDSPE PGRLEDALMH

image



Background

Alternative Names SH3 domain-binding protein 2 (3BP-2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SH3BP2 Protein
Catalog #: PEH8289

