


Product Name: Recombinant Human SMC02 Protein
Catalog #: PEH8291



Summary

Name	SMC02_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NFE2
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KMSLQMKMDC QEQLTKKNN GFFQKLVTE GAMQDLLKEI IKVDHILDRS DDEDDISSEN PQTDFLHKGM LELEAHDQD LSKQDKQETD VDEDQASTS LQFSKKNLLE LCLKGMFLKL NYWNTKIGLQ VKELGADYID GTEKIDNIIK KINVTENTVK SLLKDMLTLK GQIEKLEDRG LDLDQGTSTE VNTCNEVYEL KKKVIERLED LCKNVLLSA KLRMYQMEAE DTDSHSSEEI DTEEMEALLP QAPASFLVQK SPPRNTAWKR

image



Background

Alternative Names Single-pass membrane and coiled-coil domain-containing protein 2

Product Name: Recombinant Human SMC02 Protein
Catalog #: PEH8291



Note

For Research Use Only , Not for Diagnostic Use