


Product Name: Recombinant Human SERF1A Protein
Catalog #: PEH8292



Summary

Name	SERF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75920
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQKNMCKTQE ISKGKRKEDS LTASQRKQSS GGQKSESKMS AGPHLPLKAP RENPCFPLPA AGGSRYLAY

image



Background

Alternative Names Small EDRK-rich factor 1 (Protein 4F5) (h4F5) (SMA modifier 1)

Note

For Research Use Only , Not for Diagnostic Use