


Product Name: Recombinant Human SLC8A1 Protein
Catalog #: PEH8296



Summary

Name	NAC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32418
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRGMIIIEHEG DRPSSKTEIE MDGKVVNSHV ENFLDGALVL EVDERDQDDE EARREMARIL KELKQKHPDK EIEQLIELAN YQVLSQQQKS RAFYRIQATR LMTGAGNILK RHAADQARKA VSMHEVNTEV TENDPVSKIF FEQGTYYQCLE NCGTVALTII RRGDLTNTV FVDFRTEDGT ANAGSDYEFT EGTVVFKPGD TQKEIRVGII DDDIFEEDEN FLVHLSNVKV SSEASEDGIL

image



Background

Alternative Names Sodium/calcium exchanger 1 (Na⁺)/Ca²⁺-exchange protein 1) (Solute carrier family 8 member 1)

Note

Product Name: Recombinant Human SLC8A1 Protein
Catalog #: PEH8296



For Research Use Only , Not for Diagnostic Use