


Product Name: Recombinant Human SF3A1 Protein
Catalog #: PEH8299



Summary

Name	SF3A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15459
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVSPITGEKI PASKMQEHMR IGLDPRWLE QRDRSIREKQ SDDEVYAPGL DIESSLKQLA ERRTDIFGVE ETAIGKKIGE EEIQKPEEKV TWDGHSGSMA RTQQAAQANI TLQEQIEAIH KAKGLVPEDD TKEKIGPSKP NEIPQQPPPP SSATNIPSSA

image



Background

Alternative Names	Splicing factor 3A subunit 1 (SF3a120) (Spliceosome-associated protein 114) (SAP 114)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SF3A1 Protein
Catalog #: PEH8299

