


Product Name: Recombinant Human SF3B1 Protein
Catalog #: PEH8300



Summary

Name	SF3B1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75533
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IEAQIREIQG KKAALDEAQG VGLDSTGYDD QEIYGGSDSR FAGYVTSIAA TELEDDDDY SSSTSLLGQK KPGYHAPVAL LNDIPQSTEQ YDPFAEHRPP KIADREDEYK KHRRTMIISP ERLDPFADGG KTPDPKMNAR TYMDVMREQH LTKEEREIRQ QLAEKAKAGE LKVVNGAAAS QPPSKRKRWW DQTADQTPGA TPKKLSSWDQ AETPGHTPSL RWEDETPGRAK GSETPGATPG SKIWDPTPSH TPAGAATPGR GDTPGHATPG HGGATSSARK NRWDETPKTE

image



SDS

Background

Alternative Names	Splicing factor 3B subunit 1 (Pre-mRNA-splicing factor SF3b 155 kDa subunit) (SF3b155) (Spliceosome-associated protein 155) (SAP 155)
--------------------------	---

Product Name: Recombinant Human SF3B1 Protein
Catalog #: PEH8300



Note

For Research Use Only , Not for Diagnostic Use