


Summary

Name	ERR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11474
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSPKGSSETE TEPPVALAPG PAPTRCLPGH KEEEDGEAGAG PGEQGGGKLV LSSLPKRLCL VCGDVASGYH YGVASCEACK AFFKRTIQGS IEYSCPASNE CEITKRRRKA CQACRFTKCL RVGMLKEGVR LDRVRRGGRQK YKRRPEVDPL PFPGPFPAGP LAVAGGPRKT AAPVNALVSH LLVVEPEKLY AMPDPAGPDG HLPVATLCD LFDREIVVTI SWAKSIPGFS SLSLSDQMSV LQSVWMEVLV LGVAQRSLPL QDELAFADL VLDEEGARAA GLGELGAALL QLVRRQLQALR LEREEVLLK ALALANSDSV HIEDAEVEQ LREALHEALL EYEAGRAGPG GGAERRRAGR LLLTLPLLRQ TAGKVLAHFY GVKLEGKVPM HKLFLEMLEA

image



Background

Alternative Names Steroid hormone receptor ERR1 (Estrogen receptor-like 1) (Estrogen-related

Product Name: Recombinant Human ESRRRA Protein
Catalog #: PEH8303



receptor alpha) (ERR-alpha) (Nuclear receptor subfamily 3 group B member
1)

Note

For Research Use Only , Not for Diagnostic Use