


Product Name: Recombinant Human STX1B Protein
Catalog #: PEH8307



Summary

Name	STX1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P61266
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKDSDDEEEV VHVDRDHFMD EFFEQVEEIR GCIEKLSADV EQVKKQHSAI LAAPNPDEKT KQELEDLTAD IKKTANKVRS KLKAIEQSIE QEEGLNRSSA DLRIRKTQHS TLSRKFVEVM TEYNATQSKY RDRCKDRIQR QLEITGRITTT NEELEDMLES GKLAIFTDDI KMDSQMTKQA LNEIETRHNE IIKLETSIRE LHDMFVDMAM LVESQGEMID RIEYNVEHSV DYVERAVSDT KKAVKYQSKA

image



Background

Alternative Names Syntaxin-1B (Syntaxin-1B1) (Syntaxin-1B2)

Note

Product Name: Recombinant Human STX1B Protein
Catalog #: PEH8307



For Research Use Only , Not for Diagnostic Use