


Product Name: Recombinant Human STX3 Protein
Catalog #: PEH8308



Summary

Name	STX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13277
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQLTQDDDDTD AVEIAIDNTA FMDEFFSEIE ETRLNIDKIS EHVEEAKKLY SIILSAPIPE PKTKDDLEQL TTEIKKRANN VRNKLKSM EK HIEEDEV RSS ADLRIRKSQH SVLSRKFVEV MTKYNEAQVD FRERSKGRIQ RQLEITGKKT TDEELEEMLE SGNPAIFTSG IIDSQISKQA LSEIEGRHKD IVRLESSIKE LHDMFMDIAM LVENQGEMLD NIELNVMHTV DHVEKARDET KKAVKYQSQA RKKLIIIVL

image



Background

Alternative Names Syntaxin-3

Product Name: Recombinant Human STX3 Protein
Catalog #: PEH8308



Note

For Research Use Only , Not for Diagnostic Use