


Product Name: Recombinant Human ALYREF Protein
Catalog #: PEH8313



Summary

Name	THOC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q86V81
Host	E.coli
Species	Human
Molecular Mass	68 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DIIKLRSQR GGRGGGRGRG RAGSQGGRGG GAQAAARVNR GGGPIRNRPA IARGAAGGGG RNRPAPYSRP KQLPDKWQHD LFDSGFGGGA GVETGGKLLV SNLDFGVSDA DIQELFAEFG TLKKAHVHYD RSGRSLGTAD VHFERKADAL KAMKQYNGVP LDGRPMNIQL VTSQIDAQRR PAQSVNRGGM TRNRGAGGFG GGGGTRRGTR GGARGRGRGA GRNSKQQLSA EELDAQLDAY

image



Background

Alternative Names THO complex subunit 4 (Tho4) (Ally of AML-1 and LEF-1) (Aly/REF export factor) (Transcriptional coactivator Aly/REF) (bZIP-enhancing factor BEF)

Note

Product Name: Recombinant Human ALYREF Protein
Catalog #: PEH8313



For Research Use Only , Not for Diagnostic Use