


Product Name: Recombinant Human TOX3 Protein
Catalog #: PEH8314



Summary

Name	TOX3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15405
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDPASLDFAQ CLGYYGYSKF GNNNNYMNMA EANNAFFAAS EQTFHTPSLG DEEFEIPPIT PPPESDPALG MPDVLLPFQA LSDPLPSQGS EFTPQFPPQS LDLPSITISR NLVEQDGV LH SSSLHMDQSH TQVSQYRQDP SLIMRSIVHM TDAARSGVMP PAQLTTINQS QLSAQLGLNL GGASMPHTSP SPPASKSATP SPSSSINEED ADEANRAIGE KRAAPDSGKK PKTPKKKKKK DPNEPQKPVS AYALFFRDTQ AAIKGQNPNA TFGEVSKIVA SMWDSLGEEQ KQVYKRKTEA AKKEYLKALA

image



Background

Alternative Names TOX high mobility group box family member 3 (CAG trinucleotide repeat-containing gene F9 protein) (Trinucleotide repeat-containing gene 9 protein)

Product Name: Recombinant Human TOX3 Protein
Catalog #: PEH8314



Note

For Research Use Only , Not for Diagnostic Use