


Product Name: Recombinant Human TCEA2 Protein
Catalog #: PEH8316



Summary

Name	TCEA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15560
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IARRLDKMVT KKSAEGAMD LRELKAMPIT LHLQSTRVG MSVNALRKQS SDEEVIALAK SLIKSWKLL DASDAKARER GRGMPLPTSS RDASEAPDPS RKRPELPRAP STPRITTFPP VPVTCDAVRN KCREMLTAAL QTDHDHVAIG ADCERLSAQI EECIFRDVGN TDMKYKNRVR SRISNLKDAK NPDLRRNVLC GAITPQQIAV MTSEEMASDE LKEIRKAMTK EAIREHQMAR TGGTQTDLFT CGKCRKKNCT

image



Background

Alternative Names Transcription elongation factor A protein 2 (Testis-specific S-II) (Transcription elongation factor S-II protein 2) (Transcription elongation factor TFIIIS.I)

Product Name: Recombinant Human TCEA2 Protein
Catalog #: PEH8316



Note

For Research Use Only , Not for Diagnostic Use