

Product Name: Recombinant Human BDP1 Protein
Catalog #: PEH8317



Summary

Name	BDP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6H8Y1
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNASKATELE NKNLGPVTTA ENKDQSKLAC VHGIKGTSSIS SEVNLTERNE NQEESSQEVH MLSVAPVASS ETGPCTLGLD RGLGENSVEE PQIKDSKGDS VLTLPVPEYT PTSIPEVQQE NIINPQDLTV NLVANVPQDG EDEQAFILTL VEIPANAVEE

image



Background

Alternative Names Transcription factor TFIIIB component B'' homolog (Transcription factor IIIB 150) (TFIIIB150) (Transcription factor-like nuclear regulator)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BDP1 Protein
Catalog #: PEH8317

