

Product Name: Recombinant Human PRRG1 Protein
Catalog #: PEH8319



Summary

Name	TMG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14668
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRQTVTEGHI PFPQHLNIIT PPPPPDEVFD SGLSPGFLG YVVGSRSDSVS TRLSNCDPPP TYEATGQVN LQRSETEPHL DPPPEYEDIV

image



Background

Alternative Names Transmembrane gamma-carboxyglutamic acid protein 1 (Proline-rich gamma-carboxyglutamic acid protein 1) (Proline-rich Gla protein 1)

Note

For Research Use Only , Not for Diagnostic Use