



Product Name: Recombinant Human TRIO Protein
Catalog #: PEH8322



Summary

Name	TRIO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75962
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELNEEKRKSA RRKEFIMAEL IQTEKAYVRD LRECMDTYLW EMTSGVEEIP PGIVNKELII FGNMQEIYEF HNNIFLKELE KYEQLPEDVG HCFVTWADKF QMYVTYCKNK PDSTQLILEH AGSYFDEIQQ RHGLANSISS YLIKPVQRIT KYQLLLKELL TCCEEGKGEI KDGLEVMLSV PKRANDAMHL SMLEGFDENI ESQGELILQE SFQVWDPKTL IRKGRERHLF LFEMSLVFSK EVKDSSGRSK YLYKSKLFTS ELGVTEHVEG DPCKFALWVG RTPTSDNKIV LKASSIENKQ DWIKHIREVI QERTIHLKGA

image



Background

Alternative Names Triple functional domain protein (EC 2.7.11.1) (PTPRF-interacting protein)

Product Name: Recombinant Human TRIO Protein
Catalog #: PEH8322



Note

For Research Use Only , Not for Diagnostic Use