


Product Name: Recombinant Human TPH2 Protein
Catalog #: PEH8323



Summary

Name	TPH2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWU9
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGSSTLNKPN SGKNDDKGNK GSSKREAATE SGKTAVVFSL KNEVGGLVKA LRLFQEKRVN MVHIESRKSR RRSSEVEIFV DCECGKTEFN ELIQLLKFT TIVTLNPPEN IWTEEELED

image



Background

Alternative Names Tryptophan 5-hydroxylase 2 (EC 1.14.16.4) (Neuronal tryptophan hydroxylase) (Tryptophan 5-monoxygenase 2)

Note

For Research Use Only , Not for Diagnostic Use