


Product Name: Recombinant Human TP53I11 Protein
Catalog #: PEH8325



Summary

Name	P5I11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14683
Host	E.coli
Species	Human
Molecular Mass	8 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKHSQTDLVS RLKTRKILGV GGEDDDGEVH RSKISQVLGN EIKFTIREPL

image



Background

Alternative Names Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein)

Note

For Research Use Only , Not for Diagnostic Use