


Product Name: Recombinant Human MYO1C Protein
Catalog #: PEH8329



Summary

Name	MYO1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O00159
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQQKAVASEI FKGKKDNYPQ SVPRLFISTR LGTDEISPRV LQALGSEPIQ YAVPVVKYDR KGYKPRSRQL LLTPNAVIVV EDAKVKQRID YANLTGISVS SLSDSLFVLH VQRADNKQKG DVVLQSDHVI ETLTKTALSA NRVNSININQ GSITFAGGPG RDGTIDFTPG SELLITKAKN GHLAVVAPRL

image



Background

Alternative Names Unconventional myosin-Ic (Myosin I beta) (MMI-beta) (MMIb)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MYO1C Protein
Catalog #: PEH8329

