


Product Name: Recombinant Human UCK2 Protein
Catalog #: PEH8331



Summary

Name	UCK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZX2
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NHQQPNGGEP FLIGVSGGTA SGKSSVCAKI VQLLGQNEVD YRQKQVVILS QDSFYRVLTS EQKAKALKGQ FNFDPDAFD NELIKTLKE ITEGKTQVIP VYDFVSHSRK EETVTVYPAD VVLFEGILAF YSQEVRDLFQ MKLFVDTDAD TRLSRRVLRD ISERGRDLEQ ILSQYITFVK PAFEEFCLPT KKYADVIIPR GADNLVAINL IVQHIQDILN GGPSKRQTNG CLNGYTPSRK

image



Background

Alternative Names Uridine-cytidine kinase 2 (UCK 2) (EC 2.7.1.48) (Cytidine monophosphokinase 2) (Testis-specific protein TSA903) (Uridine monophosphokinase 2)

Note

Product Name: Recombinant Human UCK2 Protein
Catalog #: PEH8331



For Research Use Only , Not for Diagnostic Use