


**Product Name: Recombinant Human ACADVL Protein**  
**Catalog #: PEH8337**



## Summary

<b>Name</b>	ACADV_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P49748
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	29 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KSDSHPSDAL TRKKPAKAES KSFVGMFKG QLTTDQVPY PSVLNEEQTQ FLKELVEPVS RFFEEVNDPA KNDALEMVEE TTWQGLKELG AFGLQVPSEL GGVGLCNTQY ARLVEIVGMH DLGVGITLGA HQSIGFKGIL LFGTKAQKEK YLPKLASGET VAAFCLTEPS SGSDAASIRT SAVPSPCGKY YTLNGSKLWI SNGGLADIFT VFAKTPVTDP ATGAVKEKIT AFVVERGFGG ITHGPPEKMM GIKASNTAEV

## image



## Background

<b>Alternative Names</b>	Very long-chain specific acyl-CoA dehydrogenase, mitochondrial (VLCAD) (EC 1.3.8.9)
--------------------------	---

**Product Name: Recombinant Human ACADV L Protein**  
**Catalog #: PEH8337**



---

**Note**

For Research Use Only , Not for Diagnostic Use