


Product Name: Recombinant Human SORBS3 Protein
Catalog #: PEH8340



Summary

Name	VINEX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	O60504
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SEELPRSTFN YRPGAFSTVL QPSNQVLRRR EKVDNVWTEE SWNQFLQELE TGQRPKKPLV DDPGEKPSQP IEVLLERELA ELSAELDKDL RAIETRLPSP KSSPAPRRAP EQRPPAGPAS AWSSSYPHAP YLGSARSLSP HKMADGGSPF LGR RDFVYPS STRDPSASNG GGSPARREEK KRKAARLKFD FQAQSPKELT LQKGDIVYIH KEVDKNWLEG EHHGRLGIFP ANYVEVLPAD

image



Background

Alternative Names Vinexin (SH3-containing adapter molecule 1) (SCAM-1) (Sorbin and SH3 domain-containing protein 3)

Note

Product Name: Recombinant Human SORBS3 Protein
Catalog #: PEH8340



For Research Use Only , Not for Diagnostic Use