


Product Name: Recombinant Human KCNAB2 Protein
Catalog #: PEH8342



Summary

Name	KCAB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q13303
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARLSLRQTGS PGMIYSTRYG SPKRQLQFYR NLGKSGLRVS CLGLGTWVTF GGQITDEMAE QLMTLAYDNG INLFDTAEVY AAGKAEVVLG NIIKKKGWRR SSLVITTKIF WGGKAETERG LSRKHIEGL KASLERLQLE YVDVVFANRP DPNTPMEETV RAMTHVINQG MAMYWGTSRW SSMEIMEAYS VARQFNLTTP ICEQAEYHMF QREKVEVQLP ELFKIGVGA MTWSPLACGI VSGKYDSGIP PYSRASLKG Y QWLKDKILSE EGRRQQAKLK

image



Background

Alternative Names	Voltage-gated potassium channel subunit beta-2 (EC 1.1.1.-) (K(+)) channel subunit beta-2) (Kv-beta-2) (hKvbeta2)
--------------------------	---

Product Name: Recombinant Human KCNAB2 Protein
Catalog #: PEH8342



Note

For Research Use Only , Not for Diagnostic Use