


Product Name: Recombinant Human ZNF23 Protein
Catalog #: PEH8345



Summary

Name	ZNF23_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17027
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGKENVSFEL QRDFSQETDF SEASLLEKQQ EVHSAGNIKK EKSNTIDGTV KDETSPVEEC FFSQSSNSYQ CHTITGEQPS GCTGLGKSIS FDTKLVKHEI INSEERPFKC EELVEPFRCD SQLIQHQENN TEEKPYQCSE CGKAFSINEK

image



Background

Alternative Names	Zinc finger protein 23 (Zinc finger protein 359) (Zinc finger protein 612) (Zinc finger protein KOX16)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use