


Product Name: Recombinant Human ZNF354C Protein
Catalog #: PEH8346



Summary

Name	Z354C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86Y25
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PVTFRDVAVF FSQDEWLHLD SAQRALYREV MLENYSSLVS LGIPFSMPKL IHQLQQGEDP CMVEREVPSD TRLGFKTWLE TEALPHRQDI FIEETSQGMV KKESIKDGHW DINFEEAVEF ESEIEEEQEK KPLRQMIDSH EKTISEDGNH TSLELGKSLF TNTALVTQQS

image



Background

Alternative Names	Zinc finger protein 354C (Kidney, ischemia, and developmentally-regulated protein 3) (hKID3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZNF354C Protein
Catalog #: PEH8346

