


Product Name: Recombinant Human CIT Protein
Catalog #: PEH8351



Summary

Name	CTRO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14578
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQDKCHKMEQ EMTRLHRRVS EVEAVLSQKE VELKASETQR SLLEQDLATY ITECSSLKRS LEQARMEVSQ EDDKALQLLH DIREQSRKLQ EIKEQEYQAA VEEMRLMMNQ LEEDLVSARR RSDLYESELK ESRLAAEEFK RKATECQHKL LKAKDQGKPE VGEYAKLEKI NAEQQLKIQE LQEKLEKAVK ASTEATELLQ NIRQAKERA E RELEKLQNRE DSSEGIRKKL VEAEELKEKH REAQVSAQHL EVHLKQKEQH YEEKIKVLDN QIKKDLADKE TLENMMQRHE EEAHEKGKIL SEQKAMINAM DSKIRSLEQR IVELSEANKL AANSSLFTQR NMKAQEEMIS ELRQQKFYLE TQAGKLEAQN RKLLEQLEKI SHQDHSDKNR

image



Background

Alternative Names	Citron Rho-interacting kinase (CRIK) (EC 2.7.11.1) (Serine/threonine-protein
--------------------------	--

Product Name: Recombinant Human CIT Protein
Catalog #: PEH8351



kinase 21)

Note

For Research Use Only , Not for Diagnostic Use