



Product Name: Recombinant Human FANCG Protein
Catalog #: PEH8354



Summary

Name	FANCG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15287
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FAYRQGLQEL ITGNPDKALS SLHEAASGLC PRPVLVQVYT ALGSCHRKMKG NPQRALLYLV AALKEGSAWG PPLLEASRLY QQLGDTTAELE SLELLVEAL NVPCSSKAPQ FLIEVELLLP PPDLASPLHC GTQSQTKHIL ASRCLQTGRA GDAAEHYLDL LALLDSSEP RFSPPPSPPG PCMPEVFLEA AVALIQAGRA QDALTLCCEL LSRTSSLLPK MSRLWEDARK GTKELPYCPL WVSATHLLQG QAWVQLGAQK VAISEFSRCL ELLFRATPEE KEQGAAFNCE QGCKSDAALQ QLRAAALISR GLEWVASGQD TKALQDFLLS VQMCPGNRDT YFHLLQTLKR LDRRDEATAL WWRLEAQTKG SHEDALWSLP LYLESYLSWI RPSDRDAFLE EFRTSLPKSC

image



Background

Product Name: Recombinant Human FANCG Protein
Catalog #: PEH8354



Alternative Names

Fanconi anemia group G protein (Protein FCG) (DNA repair protein XRCC9)

Note

For Research Use Only , Not for Diagnostic Use