


Product Name: Recombinant Human TRIM37 Protein
Catalog #: PEH8362



Summary

Name	TRI37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94972
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSRDQHWYIT QLEAAQTSYI QQINNLKERL TIELSRTQKS RDLSPPDNHL SPQNDDALET RAKKSACSDM LLEGGPTTAS VREAKEDEED EEKIQNEDYH HELSDGDL DLVYEDEVNQ LDGSSSSASS TATSNTEEND IDEETMSGEN DVEYNNMELE EGELMEDAAA AGPAGSSHGY VGSSSRISRR THLCSAATSS LLDIDPLILI HLLDLKDRSS IENLWGLQPR PPASLLQPTA SYSRKDKDQR

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRIM37 (EC 2.3.2.27) (Mulibrey nanism protein) (RING-type E3 ubiquitin transferase TRIM37) (Tripartite motif-containing protein 37)
--------------------------	---

Product Name: Recombinant Human TRIM37 Protein
Catalog #: PEH8362



Note

For Research Use Only , Not for Diagnostic Use