


Product Name: Recombinant Human ELP1 Protein
Catalog #: PEH8363



Summary

Name	ELP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95163
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YNEALKLYSP SSQYQDISI AYGEHLMQEH MYEPAGLMFA RCGAHEKALS AFLT CGNWKQ ALCVAAQLNF TKDQLVGLGR TLAGKLVEQR KHIDAAMVLE ECAQDYEEAV LLLLEGAWE EALRLVYKYN RLDIITNVK PSILEAQKNY MAFLDSQTAT FSRHKKRLLV VRELKEQAQQ AGLDDEVPHG QESDLFSETS SVVSGSEMSG KYSHSNSRIS ARSSKNRRKA ERKKHSLKEG SPLEDLALLE ALSEVVQNTE NLKDEVYHIL KVLFLFEFDE QGRELQKAFE DTLQLMERSL PEIWTLYQQ NSATPVLGPN STANSIMASY

image



Background

Alternative Names Elongator complex protein 1 (ELP1) (IkappaB kinase complex-associated protein) (IKK complex-associated protein) (p150)

Product Name: Recombinant Human ELP1 Protein
Catalog #: PEH8363



Note

For Research Use Only , Not for Diagnostic Use