


Product Name: Recombinant Human RPL13A Protein
Catalog #: PEH8370



Summary

Name	RL13A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P40429
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAALKVVRLK PTRKFAYLGR LAHEVGWKYQ AVTATLEEKR KEKAKIHYRK KKQLMRLRKQ AEKNVEKKID KYTEVLKTHG

image



Background

Alternative Names	60S ribosomal protein L13a (23 kDa highly basic protein) (Large ribosomal subunit protein uL13)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use