


Product Name: Recombinant Human GSTM5 Protein
Catalog #: PEH8372



Summary

Name	GSTM5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P46439
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGLAHAIRLL LEYTDSSYVE KKYTLGDAPD YDRSQWLNEK FKLGLDFPNL PYLIDGAHKI TQSNAILRYI ARKHNLCGET EEEKIRVDIL ENQVMDNHME LVRLCYDPDF EKLKPKYLEE LPEKLYSE FLGKRPWFAG DKITFVDFLA YDVLDMKRIF EPKCLDAFLN LKDFISRFEG LKKISAYMKS SQFLRGLLFG

image



Background

Alternative Names Glutathione S-transferase Mu 5 (EC 2.5.1.18) (GST class-mu 5) (GSTM5-5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GSTM5 Protein
Catalog #: PEH8372

