


Product Name: Recombinant Human CENPF Protein
Catalog #: PEH8373



Summary

Name	CENPF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49454
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLPTRALQKI QELEGQLDKL KKEKQQRQFQ LDSLEAALQK QKQKVENEKT EGTNLKRENQ RLMEICESLE KTKQKISHL QVKESQVNFQ EGQLNSGKKQ IEKLEQELKR CKSELEERSQQ AAQSADVSLN PCNTPQKIFT TPLTPSQYYS GSKYEDLKEK YNKEVEERKR LEAEVKALQA KKASQTLPPA TMNHRDIARH QASSSVFSWQ QEKTPSHLSS NSQRTPIRRD FSASYFSGEQ EVTPSRSTLQ IGKRDANSSF FDNSSSPHLL DQLKAQNQEL RNKINELELR LQGHEKEMKG

image



Background

Alternative Names	Centromere protein F (CENP-F) (AH antigen) (Kinetochore protein CENPF) (Mitotin)
--------------------------	--

Product Name: Recombinant Human CENPF Protein
Catalog #: PEH8373



Note

For Research Use Only , Not for Diagnostic Use