


Product Name: Recombinant Human TTPA Protein
Catalog #: PEH8374



Summary

Name	TTPA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P49638
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPQLNALPDH SPLLQPLGAA LRRRAREAGV PLAPLPLTDS FLLRFLRARD FDLDLAWRL KNYWKWRAEC PEISADLHPR SIIGLLKAGY HGVLRSDPT GSKVLIYRIA HWDPKVFTAY DVFRVSLITS ELIVQEVETQ RNGIKAIFDL EGWQFSHAFQ ITPSVAKKIA AVLTDSFPLK VRGIHLINEP VIFHAVFSMI KPFLTEKIKE RIHMHGNNYK QSLLQHFPDI LPLEYGGEEF

image



Background

Alternative Names Alpha-tocopherol transfer protein (Alpha-TTP)

Note

Product Name: Recombinant Human TTPA Protein
Catalog #: PEH8374



For Research Use Only , Not for Diagnostic Use