


Product Name: Recombinant Human BRCA2 Protein
Catalog #: PEH8376



Summary

| | |
|------------------------------|--|
| Name | BRCA2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | P51587 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 36 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | QNISEKDLLD TENKRKKDFL TSENSLPRIS SLPKSEKPLN EETVVNKRDE EQHLESHTDC ILAVKQAISG TSPVASSFQG IKKSIFRIRE SPKETFNASF SGHMTDPNFK KETEASESGL EIHTVCSQKE DSLCPNLIDN GSWPATTQTN SVALKNAGLI STLKKTNKF IYAIHDETSY KGKKIPKDQK SELINCSAQF EANAPEAPLT FANADSGLLH SSVKRSCSQN DSEPTLSLT SSFGTILRKC SRNETCSNNT VISQDLDYKE AKCNKEKLQL FITPEADSL SCLQEGQCEND PKSKKVSDIK |

image



Background

Alternative Names Breast cancer type 2 susceptibility protein (Fanconi anemia group D1 protein)

Product Name: Recombinant Human BRCA2 Protein
Catalog #: PEH8376



Note

For Research Use Only , Not for Diagnostic Use