


Product Name: Recombinant Human GEMIN4 Protein
Catalog #: PEH8378



Summary

Name	GEMI4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57678
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFQEDLNTTF NQLTQSASEQ GLAKAVASVA RLVIVHPEVT VKKMCSLAVV NLGTHKFLAQ ILTAFPALRF VEEQGPNSSA TFMVSLKET VWMKFSTPKE EKQFLELLNC LMSPVKPQGI PVAALLEPDE VLKEFVLPFL RLDVEEVDLS LRIFIQTLEA

image



Background

Alternative Names Gem-associated protein 4 (Gemin-4) (Component of gems 4) (p97)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GEMIN4 Protein
Catalog #: PEH8378

