


Summary

Name	RBX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P62877
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGTNSGAGKK RFEVKKWNAV ALWAWDIVVD NCAICRNHIM DLCIECQANQ ASATSEECTV AWGVCNHAFH FHCISRWLKT RQVCPLDNRE

image



Background

Alternative Names	E3 ubiquitin-protein ligase RBX1 (EC 2.3.2.27) (EC 2.3.2.32) (E3 ubiquitin-protein transferase RBX1) (Protein ZYP) (RING finger protein 75) (RING-box protein 1) (Rbx1) (Regulator of cullins 1) (ROC1) [Cleaved into: E3 ubiquitin-protein ligase RBX1, N-terminally processed (E3 ubiquitin-protein transferase RBX1, N-terminally processed)]
--------------------------	--

Note

Product Name: Recombinant Human RBX1 Protein
Catalog #: PEH8379



For Research Use Only , Not for Diagnostic Use