


Product Name: Recombinant Human RASSF7 Protein
Catalog #: PEH8381



Summary

Name	RASF7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02833
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RVQRNAEELG HEAFWEQELR REQAREREGQ ARLQALSAAT AEHAARLQAL DAQARALEAE LQLAAEAPGP PSPMASATER LHQDLAVQER QSAEVQGSLA LVSRALEAAE RALQAQAQEL EELNRELRC NLQQFIQQTG AALPPPPRPD RGPPGTQGPL PPAREESLLG APSESHAGAQ PRPRGGPHDA

image



Background

Alternative Names	Ras association domain-containing protein 7 (HRAS1-related cluster protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RASSF7 Protein
Catalog #: PEH8381

