


Product Name: Recombinant Human KIF5A Protein
Catalog #: PEH8383



Summary

Name	KIF5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12840
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QWKKKYEKEK EKTKAQKETI AKLEAELSRW RGENVPETE RLAGEEAALG AELCEETPVN DNSSIVVRIA PEERQKYEEE IRRLYQLDD KDDEINQQSQ LIEKCLKQML DQEELLVSTR GDNEKVQREL SHLQSENDAA KDEVKEVLQA LEELAVNYDQ KSQEVVEKSQ

image



Background

Alternative Names Kinesin heavy chain isoform 5A (Kinesin heavy chain neuron-specific 1) (Neuronal kinesin heavy chain) (NKHC)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KIF5A Protein
Catalog #: PEH8383

