


Summary

Name	GNDS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12967
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LINGVIYSIS LRKVQLHHGG NKGQRWLGYE NESALNLYET CKVRTVKAGT LEKLVEHLVP AFQGSLSYV TIFLCTYRAF TTTQQVLDLL FKRYGRCDAL TASSRYGCIL PYSDEDGGPQ DQLKNAISSI LGTWLDQYSE DFCQPPDFPC LKQLVAYVQL NMPGSDLERR AHLLLAQLEH SEPIEAEPEA LSPVPALKPT

image



Background

Alternative Names	Ral guanine nucleotide dissociation stimulator (RalGDS) (Ral guanine nucleotide exchange factor) (RalGEF)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RALGDS Protein
Catalog #: PEH8384

