


**Product Name: Recombinant Human REST Protein**  
**Catalog #: PEH8385**



## Summary

<b>Name</b>	REST_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q13127
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	36 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	PFKCDQCSYV ASNQHEVTRH ARQVHNGPKP LNCPHCDYKT ADRSNFKKHV ELHVNPRQFN CPVCDYAASK KCNLQYHFKS KHPTCPNKTM DVSKVKLKKT KKREADLPDN ITNEKTEIEQ TKIKGDVAGK KNEKSVKAEK RDVSKEKKPS NNVSVIQVTT RTRKSVTEVK EMDVHTGSNS EKFSKTKKSK RKLEVDSHSL HGPVNDEESS TTKKKKVESK SKNNSQEVPK GDSKVEENKK QNTCMKKSTK KKTLKNKSSK KSSKPPQKEP VEKGSAQMDP PQMGAPTEA VQKGPVQVEP PPPMEHAQME

## image



## Background

<b>Alternative Names</b>	RE1-silencing transcription factor (Neural-restrictive silencer factor) (X2 box repressor)
--------------------------	--

**Product Name: Recombinant Human REST Protein**  
**Catalog #: PEH8385**

---



**Note**

For Research Use Only , Not for Diagnostic Use