

Product Name: Recombinant Human SCRIB Protein
Catalog #: PEH8388



Summary

Name	SCRIB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14160
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEEEEGSPQ EEEEEEEEN RAEEEEASTE EEDKEGAVVS APSVKGVSFD QANNLLIEPA RIEEEELTLT ILRQTGGLGI SIAGGKGSTP YKGDDEGIFI SRVSEEGPAA RAGVRVGDKL LEVNGVALQG AEHHEAVEAL RGAGTAVQMR VWRERMVEPE NAVTITPLRP EDDYSPRERR GGGLRLPLL PESPGPLRQR HVACLARSER GLGFSIAGGK GSTPYRAGDA GIFVSRIAEG GAAHRAGTLQ VGDRVLSING VDVTEARHDH AVSLLTAASP TIALLLEREA GGPLPPSPLP HSSPPTAAVA TTSITTATPG VPGLPSLAPS LLAAALEGPY PVEEIRLPRA GGPLGLSIVG GSDHSSHPFG VQEPGVFISK VLPRGLAARS GLRVGDRILA VNGQDVRDAT HQEAVSALLR PCLELSLLVR RDPAPPGLRE LCIQKAPGER LGISIRGGAR GHAGNPRDPT DEGIFISKVS PTGAAGRDGR LRVGLRLLEV NQQSLLGLTH GEAVQLLRSV GDTLTLVLCD

image



Background

Product Name: Recombinant Human SCRIB Protein
Catalog #: PEH8388



Alternative Names

Protein scribble homolog (Scribble) (hScrib) (Protein LAP4)

Note

For Research Use Only , Not for Diagnostic Use