



Product Name: Recombinant Human HLTf Protein
Catalog #: PEH8389



Summary

Name	HLTF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14527
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTEFDKLFED LKEDDKTHEM EPAEAIETPL LPHQKQALAW MVSRENSKEL PPFWEQRNDL YYNTITNFSE KDRPENVHGG ILADDMGLGK TLTAIAVILT NFHDGRPLPI ERVKKNLLKK EYNVNDDSMK LGGNNTSEKA

image



Background

Alternative Names

Helicase-like transcription factor (EC 2.3.2.27) (EC 3.6.4.-) (DNA-binding protein/plasminogen activator inhibitor 1 regulator) (HIP116) (RING finger protein 80) (RING-type E3 ubiquitin transferase HLTf) (SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 3) (Sucrose nonfermenting protein 2-like 3)

Product Name: Recombinant Human HLTF Protein
Catalog #: PEH8389



Note

For Research Use Only , Not for Diagnostic Use