


Product Name: Recombinant Human CHD4 Protein
Catalog #: PEH8390



Summary

Name	CHD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14839
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EMGEEEEVER EIIKQESVD PDYWEKLLRH HYEQQQEDLA RNLGKGRIR KQVNYNDGSQ EDRDWQDDQS DNQSDYSVAS EEGDEDFDER SEAPRRPSRK GLRNDKDKPL PLLARVGGN IEVLGFNARQ RKAFLNAIMR YGMPPQDAFT TQWLVRDLRG KSEKEFKAYV SLFMRHLCEP GADGAETFAD GVPREGLSRQ HVLTRIGVMS LIRKKVQEFE HVNGRWSMPE LAEVEENKKM SQPGSPSPKT PTPSTPGDTQ

image



Background

Alternative Names Chromodomain-helicase-DNA-binding protein 4 (CHD-4) (EC 3.6.4.12) (ATP-dependent helicase CHD4) (Mi-2 autoantigen 218 kDa protein) (Mi2-beta)

Product Name: Recombinant Human CHD4 Protein
Catalog #: PEH8390



Note

For Research Use Only , Not for Diagnostic Use