


Product Name: Recombinant Human SIX1 Protein
Catalog #: PEH8393



Summary

Name	SIX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15475
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PREKRELAEA TGLTTTQVSN WFKNRRQRDR AAERAKERENT ENNNSSSNKQ NQLSPLEGGK PLMSSSEEEF SPPQSPDQNS VLLQGNMGH ARSSNYSLPG LTASQPSHGL QTHQHQLQDS LLGPLTSSLV

image



Background

Alternative Names Homeobox protein SIX1 (Sine oculis homeobox homolog 1)

Note

For Research Use Only , Not for Diagnostic Use