


Product Name: Recombinant Human TSC22D1 Protein
Catalog #: PEH8394



Summary

Name	T22D1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15714
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVNTSRFRVV KLDSSSEPFK KGRWTCTEFY EKENAVPATE GVLINKVVET VKQNPIEVTS ERESTSGSSV SSSVSTLSHY TESVGSSEMGA APTVVVQQQQ QQQQQQQQQP ALQGVTLQQM DFGSTGPQSI PAVSIPQSI QSQISQVQLQ SQELSYQQKQ GLQPVPLQAT MSAATGIQPS PVNVVGVVTS LGQQPSISL AQPQLPYSQA APPVQTPLPG APPPQQQLQYG QQQPMVSTQM APGHVKSVTQ NPASEYVQQQ PILQTAMSSG QPSSAGVGAG TTVIPVAQPQ GIQLPVQPTA VPAQPAGASV

image



Background

Alternative Names TSC22 domain family protein 1 (Cerebral protein 2) (Regulatory protein TSC-22) (TGFB-stimulated clone 22 homolog) (Transforming growth factor

Product Name: Recombinant Human TSC22D1 Protein
Catalog #: PEH8394



beta-1-induced transcript 4 protein)

Note

For Research Use Only , Not for Diagnostic Use