


Product Name: Recombinant Human KIF7 Protein
Catalog #: PEH8398



Summary

Name	KIF7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2M1P5
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEEIAAFQRK RRSNGSVV SLEQQKIEE QKKWLDQEME KVLQRRRALE ELGEELHKRE AILAKKEALM QEKTGLESKR LRSSQALNED IVRVSSRLEH LEKELSEKSG QLRQGSASQ QQIRGEIDSL RQEKDSLKKQ RLEIDGKLRQ GSLSPPEER TLFQLDEAIE ALDAAIEYKN EAITCRQRL

image



Background

Alternative Names Kinesin-like protein KIF7

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KIF7 Protein
Catalog #: PEH8398

