



Product Name: Recombinant Human AGBL2 Protein
Catalog #: PEH8400



Summary

Name	CBPC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5U5Z8
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKFTQCLAEL KELLRQEIHK KFHELQDVD LEGSWSDISL SDIESSTSGS DSSLSDGLPV HLANIADELT QKKKMFKKKK KKS LQTRKQR NEQYQKKNLM QKLKLTEDTS EKAGFASTLQ KQPTFFKNSE NSSFLPMKNE NPRLNETNLN RRDKDTPLDP SMATLILPKN KGRMQNKKPG FTVSCSPKRT INSSQEPAPG

image



Background

Alternative Names Cytosolic carboxypeptidase 2 (EC 3.4.17.-) (ATP/GTP-binding protein-like 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AGBL2 Protein
Catalog #: PEH8400

