


Summary

Name	SDCG8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86SQ7
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NFEKTKALIQ CDQLRKELER QAERLEKELA SQQEKRAIEK DMMKKEITKE REYMGSKMLI LSQNIQLEA QVEKVTKEKI SAINQLEEQ SQLASREMDV TKVCGEMRYQ LNKTNMEKDE AEKEHREFRA KTNRDLEIKD QEIEKLRIEL DESKQHLEQE QQKAALAREE CLRLTELLGE SEHQLHLTRQ EKDSIQQSFS KEAKAALQA QQREQELTQK IQQMEAQHDK TENEQYLLLT SQNTFLTCLK EECCTLAKKL EQISQKTRSE IAQLSQEKRY TYDKLGKLQR RNEELEEQCV QHGRVHETMK QRLRQLDKHS

image



Background

Alternative Names	Serologically defined colon cancer antigen 8 (Antigen NY-CO-8) (Centrosomal colon cancer autoantigen protein) (hCCCAP)
--------------------------	--

Product Name: Recombinant Human SDCCAG8 Protein
Catalog #: PEH8403



Note

For Research Use Only , Not for Diagnostic Use