


Product Name: Recombinant Human RPTOR Protein
Catalog #: PEH8407



Summary

Name	RPTOR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N122
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSTSSSLTN DVAKQPVSRD LPSGRPGTTG PAGAQTTPHS HQFPRTKMF DKGPEQTADD ADDAAGHKSF ISATVQTGFC DWSARYFAQP VMKIPEEHD

image



Background

Alternative Names	Regulatory-associated protein of mTOR (Raptor) (p150 target of rapamycin (TOR)-scaffold protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use